

BF-2000 Solvent Evaporator

Features:

- * Handling multi samples simultaneously, saving time and labors;
- * A needle valve adjusts gas flow at each position;
- * Easy operation
- * The level of a purge plate can be easily moved up or down for different sample tubs;
- * Wide selection: 15-tubes, 30-tubes, 96-well or customized module



Specification

Catalog No.	BF2000-F	BF2000-15AK	BF2000-30AK
Description	15-Well Fix Volumn Solvent Evaporator	15-Well Solvent Evaporator	30-Wellt Solvent Evaporator
No. of Sample Wells	15	15	30
Flow Adjustable at Each	YES	YES	YES
Well	0-5L/M	0-5L/M	0-8L/M
Gas Flow Rate	310x220x420	310x220x430	310x220x430
Overall Dimentions	9kg	8kg	9kg
(LxWxH)mm	200w	560w	560w
Weight	110v 60Hz or 220v 50Hz	110v 60Hz or 220v 50Hz	110v 60Hz or 220v 50Hz
Power			
Well OD (could be customized)	φ16mm	φ16mm	φ16mm
Inlet Pressure	<60psi	<60psi	<60psi
Volumn left	0.5ml	na	na
Temperature range	ambient - 150°C ±2°C Cool: 0°C ±3°C	ambient - 150°C ±2°C	ambient - 150°C ±2°C
Unit Price (USD)	5,000.00	2,000.00	2,500.00

Evaporator Accessories

Catalog No.	Description	Unit Price (USD)
BF2000-801A	BF2000 adjustable Purge Needle 1/PK	6.00
BF2000-T0	Adjustable valve	12.00
BF2000-T2	BF2000-15AK model purge plate(not-incl. needle)	260.00
BF2000-T4	BF2000-30AK model purge plate(not-incl. needle)	360.00
BF2000-3	Standard 15 well heating block	120.00
BF2000-30	Standard 30 well heating block	200.00
BF2000-4	Custom 15 well heating block	160.00
BF2000-40	Custom 30 well heating block	280.00
BF2000-6	Flowmeter 0.5-5L/min	70.00
BF2000-96	96 well purge head	650
BF2000-7	Smart table 7411D	60.00

